


Molecular
Dimensions



PACT *premier*[™] HT-96 Eco Screen

PACT *premier* is a pH, Anion, Cation crystallization trial devised to test pH within a PEG/Ion screen environment.

MD1-36-ECO

MD1-36-ECO contains 96 x 1 mL cacodylate-free reagents in a deep-well block.

Features of PACT *premier*

- A modern, comprehensive PEG/ion screen - the most effective systematic screen available to date.
-
- This 96-well screen is really 3 screens in one:
 - 24-well pH/PEG screen
 - 24-well cation/PEG screen
 - 48-well anion/PEG screen
- Cacodylate replaced with MES.

Rationale for a new PEG/Ion screen

The first step in crystallization is often to reach for a commercially available “sparse matrix” kit, and hope that one of the conditions produces something that looks harvestable, or optimizable. If no obvious leads come out of the screen, it is hard to learn anything from the negative (precipitate and clear) results.

There are a few screens that try to test crystallization space in a more rational manner – for example, the Clear Strategy Screen and The Solubility Tool Kit.




Traditional PEG/Ion screens provide a logical test of seven cations and eleven anions using PEG 3350 as the precipitation agent. However, the user has no control over pH and hence cannot determine the effect of one cation or anion over another.

For this reason the PEG/ION/pH (PACT) screen has been developed to systematically test the effect of pH, anions and cations, using PEG as the precipitant. This screen has been implemented very successfully at the Netherlands Cancer Institute (NKI), and at the Oxford Protein Production Facility (OPPF).



PACT *premier*

pH/PEG screen

This consists of four broad range buffer systems (1) versus PEG 1500. These buffers allow one to scan the pH range from 4 to 9, without changing the chemistry of the system, so effectively isolating the effect of pH from the effect of the buffer that causes the change in pH.



Protein crystals grown successfully with PACT *premier*.



PACT premier

Cation/PEG screen

This is made up of six cations (all with chloride counter ions) that are combined with PEG 6000 at four different pHs: Acetate pH 5, MES pH 6, HEPES pH 7 and Tris pH 8. The cations tested are Na⁺, NH₄⁺, Li⁺, Mg²⁺, Ca²⁺ and Zn²⁺. The zinc ion is tested at lower concentration than the other cations in the screen (0.01 M vs. 0.2 M)

Anion/PEG screen

This is made up of 12 anions, with either sodium or potassium counter ions, which are tested at 0.2 M against PEG 3350. The anions include fluoride, bromine, iodide, thiocyanate, nitrate, formate, acetate, sulfate, tartrate, phosphate, citrate and malonate. The phosphate solution is tested at a concentration of 0.02 M. Chloride is not included here as it is the counter ion in the cation screen. Three sets of reagents are tested at pH 6.5, 7.5, and 8.5 with the Bis-Tris-Propane buffer system whilst one set of reagents is tested without buffering.

References:

(1) Newman *et al* (2005). Towards rationalization of crystallization screening for small- to medium-sized academic laboratories: the PACT/JCSG+ strategy. *Acta Cryst.* **D61**, 1426-1431.

Formulation Notes

PACT premier reagents are formulated using ultrapure water (>18.0 MΩ) and are sterile-filtered using 0.22 μm filters. No preservatives are added.

Final pH may vary from that specified on the datasheet. Molecular Dimensions will be happy to discuss the precise formulation of individual reagents.


Individual reagents and stock solutions for optimization are available from Molecular Dimensions.

Enquiries regarding **PACT-premier** formulation, interpretation of results or optimization strategies are welcome. Please e-mail, fax or phone your query to Molecular Dimensions.


Contact and product details can be found at www.moleculardimensions.com

PACT premier was developed by Janet Newman, and was tested in the laboratory of Anastassis Perrakis at the Netherlands Cancer Institute as part of the SPINE programme and is manufactured under license by Molecular Dimensions.

PACT premier when used together with **JCSG plus** as a primary screening strategy is an extremely powerful and successful combination, (i.e. a combination of a modern sparse matrix approach and an information yielding systematic trial).



Molecular
Dimensions




PACT premier HT-96 Eco Screen


Conditions A1- D12

MD1-36-ECO

Well #	Conc. Salt	Conc. Buffer	pH	Conc. Precipitant
A1		0.1 M SPG	4.0	25 % w/v PEG 1500
A2		0.1 M SPG	5.0	25 % w/v PEG 1500
A3		0.1 M SPG	6.0	25 % w/v PEG 1500
A4		0.1 M SPG	7.0	25 % w/v PEG 1500
A5		0.1 M SPG	8.0	25 % w/v PEG 1500
A6		0.1 M SPG	9.0	25 % w/v PEG 1500
A7	0.2 M Sodium chloride	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
A8	0.2 M Ammonium chloride	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
A9	0.2 M Lithium chloride	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
A10	0.2 M Magnesium chloride hexahydrate	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
A11	0.2 M Calcium chloride dihydrate	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
A12	0.01 M Zinc chloride	0.1 M Sodium acetate	5.0	20 % w/v PEG 6000
B1		0.1 M MIB	4.0	25 % w/v PEG 1500
B2		0.1 M MIB	5.0	25 % w/v PEG 1500
B3		0.1 M MIB	6.0	25 % w/v PEG 1500
B4		0.1 M MIB	7.0	25 % w/v PEG 1500
B5		0.1 M MIB	8.0	25 % w/v PEG 1500
B6		0.1 M MIB	9.0	25 % w/v PEG 1500
B7	0.2 M Sodium chloride	0.1 M MES	6.0	20 % w/v PEG 6000
B8	0.2 M Ammonium chloride	0.1 M MES	6.0	20 % w/v PEG 6000
B9	0.2 M Lithium chloride	0.1 M MES	6.0	20 % w/v PEG 6000
B10	0.2 M Magnesium chloride hexahydrate	0.1 M MES	6.0	20 % w/v PEG 6000
B11	0.2 M Calcium chloride dihydrate	0.1 M MES	6.0	20 % w/v PEG 6000
B12	0.01 M Zinc chloride	0.1 M MES	6.0	20 % w/v PEG 6000
C1		0.1 M PMTP	4.0	25 % w/v PEG 1500
C2		0.1 M PMTP	5.0	25 % w/v PEG 1500
C3		0.1 M PMTP	6.0	25 % w/v PEG 1500
C4		0.1 M PMTP	7.0	25 % w/v PEG 1500
C5		0.1 M PMTP	8.0	25 % w/v PEG 1500
C6		0.1 M PMTP	9.0	25 % w/v PEG 1500
C7	0.2 M Sodium chloride	0.1 M HEPES	7.0	20 % w/v PEG 6000
C8	0.2 M Ammonium chloride	0.1 M HEPES	7.0	20 % w/v PEG 6000
C9	0.2 M Lithium chloride	0.1 M HEPES	7.0	20 % w/v PEG 6000
C10	0.2 M Magnesium chloride hexahydrate	0.1 M HEPES	7.0	20 % w/v PEG 6000
C11	0.2 M Calcium chloride dihydrate	0.1 M HEPES	7.0	20 % w/v PEG 6000
C12	0.01 M Zinc chloride	0.1 M HEPES	7.0	20 % w/v PEG 6000
D1		0.1 M MMT	4.0	25 % w/v PEG 1500
D2		0.1 M MMT	5.0	25 % w/v PEG 1500
D3		0.1 M MMT	6.0	25 % w/v PEG 1500
D4		0.1 M MMT	7.0	25 % w/v PEG 1500
D5		0.1 M MMT	8.0	25 % w/v PEG 1500
D6		0.1 M MMT	9.0	25 % w/v PEG 1500
D7	0.2 M Sodium chloride	0.1 M Tris	8.0	20 % w/v PEG 6000
D8	0.2 M Ammonium chloride	0.1 M Tris	8.0	20 % w/v PEG 6000
D9	0.2 M Lithium chloride	0.1 M Tris	8.0	20 % w/v PEG 6000
D10	0.2 M Magnesium chloride hexahydrate	0.1 M Tris	8.0	20 % w/v PEG 6000
D11	0.2 M Calcium chloride dihydrate	0.1 M Tris	8.0	20 % w/v PEG 6000
D12	0.002 M Zinc chloride	0.1 M Tris	8.0	20 % w/v PEG 6000



Molecular
Dimensions





PACT premier HT-96 Eco Screen

Conditions E1- H12

MD1-36-ECO

Well #	Conc. Salt	Conc. Buffer	pH	Conc. Precipitant
E1	0.2 M Sodium fluoride			20 % w/v PEG 3350
E2	0.2 M Sodium bromide			20 % w/v PEG 3350
E3	0.2 M Sodium iodide			20 % w/v PEG 3350
E4	0.2 M Potassium thiocyanate			20 % w/v PEG 3350
E5	0.2 M Sodium nitrate			20 % w/v PEG 3350
E6	0.2 M Sodium formate			20 % w/v PEG 3350
E7	0.2 M Sodium acetate trihydrate			20 % w/v PEG 3350
E8	0.2 M Sodium sulfate			20 % w/v PEG 3350
E9	0.2 M Potassium sodium tartrate tetrahydrate			20 % w/v PEG 3350
E10	0.02 M Sodium/potassium phosphate			20 % w/v PEG 3350
E11	0.2 M Sodium citrate tribasic dihydrate			20 % w/v PEG 3350
E12	0.2 M Sodium malonate dibasic monohydrate			20 % w/v PEG 3350
F1	0.2 M Sodium fluoride	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F2	0.2 M Sodium bromide	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F3	0.2 M Sodium iodide	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F4	0.2 M Potassium thiocyanate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F5	0.2 M Sodium nitrate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F6	0.2 M Sodium formate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F7	0.2 M Sodium acetate trihydrate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F8	0.2 M Sodium sulfate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F9	0.2 M Potassium sodium tartrate tetrahydrate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F10	0.02 M Sodium/potassium phosphate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F11	0.2 M Sodium citrate tribasic dihydrate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
F12	0.2 M Sodium malonate dibasic monohydrate	0.1 M Bis-Tris propane	6.5	20 % w/v PEG 3350
G1	0.2 M Sodium fluoride	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G2	0.2 M Sodium bromide	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G3	0.2 M Sodium iodide	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G4	0.2 M Potassium thiocyanate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G5	0.2 M Sodium nitrate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G6	0.2 M Sodium formate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G7	0.2 M Sodium acetate trihydrate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G8	0.2 M Sodium sulfate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G9	0.2 M Potassium sodium tartrate tetrahydrate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G10	0.02 M Sodium/potassium phosphate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G11	0.2 M Sodium citrate tribasic dihydrate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
G12	0.2 M Sodium malonate dibasic monohydrate	0.1 M Bis-Tris propane	7.5	20 % w/v PEG 3350
H1	0.2 M Sodium fluoride	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H2	0.2 M Sodium bromide	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H3	0.2 M Sodium iodide	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H4	0.2 M Potassium thiocyanate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H5	0.2 M Sodium nitrate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H6	0.2 M Sodium formate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H7	0.2 M Sodium acetate trihydrate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H8	0.2 M Sodium sulfate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H9	0.2 M Potassium sodium tartrate tetrahydrate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H10	0.02 M Sodium/potassium phosphate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H11	0.2 M Sodium citrate tribasic dihydrate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350
H12	0.2 M Sodium malonate dibasic monohydrate	0.1 M Bis-Tris propane	8.5	20 % w/v PEG 3350



Abbreviations:

HEPES; N-(2-hydroxyethyl)-piperazine-N'-2-ethanesulfonic acid, **MES**; 2-(N-morpholino)ethanesulfonic acid, **PEG**; Polyethylene glycol, **Tris**; 2-Amino-2-(hydroxymethyl)propane-1,3-diol, **SPG buffer**; Succinic Acid, Phosphate, Glycine, **MIB buffer**; Malonic acid, Imidazole, Boric acid, **PMTP buffer**; Propionic acid, MES, Bis-tris propane, **MMT buffer**; Malic acid, MES, Tris.


Manufacturer's safety data sheets are available from our website or by scanning the QR code here:




Re-Ordering details:

Catalogue Description	Pack size	Catalogue Code
PACT <i>premier</i> [™]	96 x 10 mL	MD1-29
PACT <i>premier</i> [™] HT-96	96 x 1 mL	MD1-36
PACT <i>premier</i> [™] FX-96	96 x 100 µL	MD1-36-FX
Eco Screens		
PACT <i>premier</i> [™] Eco Screen	96 x 10 mL	MD1-29-ECO
PACT <i>premier</i> [™] HT-96 Eco Screen	96 x 1 mL	MD1-36-ECO
Green Screens (contain fluorescent green dye - ideal for UV)		
PACT <i>premier</i> [™] Green Screen	96 x 10 mL	MD1-55
PACT <i>premier</i> [™] HT-96 Green Screen	96 x 1 mL	MD1-52
Combo Packs		
Super2 Combo Value Pack (JCSG- <i>plus</i> [™] + PACT <i>premier</i> [™])	2 x 96 x 10 mL	MD1-75
Super2 Combo HT-96 Value Pack (JCSG- <i>plus</i> [™] HT-96 + PACT <i>premier</i> [™] HT-96)	2 x 96 x 10 mL	MD1-75-HT
Single Reagents		
PACT <i>premier</i> [™] single reagents	100 mL	MDSR-29-tube number
PACT <i>premier</i> [™] HT-96 single reagents	100 mL	MDSR-36-well number

For PACT *premier*[™] stock solutions please visit the Optimization section on our website.



Molecular Dimensions



JCSG-plus™ HT-96 Eco screen



MD1-40-ECO

JCSG-plus™ is the screen of choice for initial screening experiments.
The most complete sparse matrix screen available today.

MD1-40-ECO is presented as 96 x 1 mL cacodylate-free and dioxane-free conditions in a deep-well block.

Features of JCSG-plus™ HT-96:


- Optimized sparse matrix screen without cacodylate (MES used instead) or dioxane.
- Reduced redundancy.
- Screens classic PEG and salt conditions.
- Access more areas of crystallization space.
- Neutralised organic acids: Formate, acetate, citrate, succinate, malate, malonate.
- More organic and polyalcohol conditions
- Precipitant synergy.
- Wide pH range 4.0 – 10.0.





Introduction

Commercially available sparse matrix screens are devised using conditions based on previously successful crystallization conditions. Since increasing numbers of researchers now use commercially available sparse matrix screens, the same sub-sets of conditions are used repeatedly. A number of structural genomics initiatives have published results of data-mining exercises using internally consistent datasets and analysing negative results as well as positive hits. The results have been startling!

Members of the Joint Centre for Structural Genomics analysed the crystallization of over 500 different proteins against commercially available sparse matrix screens totalling 480 conditions, compiled to sample a wide range of precipitant, buffer, additive and pH. The **core screen (JCSG)** was developed when data mining revealed massive redundancy between clusters of conditions in commercial screens, particularly where high molecular weight PEGs are used as precipitants (1). Using a novel algorithm, members of the JCSG identified "conditions most essential for promoting crystal formation for the most diverse set of proteins. **JCSG-plus™ supersedes the JCSG Core Screen and Index screens.**



Analysis of precipitants used in JCSG-plus™



In-filling the optimized screen

The second issue to come to light was that even extensive suites of sparse matrix screens represent incomplete coverage of crystallisation space – 480 conditions failed to crystallise 15% of the target proteins.

The **JCSG-*plus***[™] screen is supplemented with additional conditions to provide a more complete coverage of crystallisation space and improved chemical complementarity (2).

- i. In-filling the pH profile
- ii. introduce conditions using neutralised organic acids as the precipitant (3)
- iii. expanded range of organic and polyalcohol conditions
- iv. precipitant synergy

Usage

JCSG-*plus*[™] is designed for the rapid, efficient screening for crystallization leads of a new protein or preparation. In the first instance, drops should be set-up using equal volumes of protein solution and reagent. Protein samples should be in a minimal solvent system containing a low concentration of buffer. Starting protein concentrations should be between 5 mg/ml and 40 mg/ml. Protein concentration can be varied in subsequent rounds depending on initial results.

The conditions in JCSG-*plus*[™] are compatible with all commonly used crystallisation methods, sitting drop, hanging drop, sandwich drop, microbatch, vapour microbatch and microdialysis.

The JCSG-*plus*[™] sparse matrix screen is highly effective when used alongside a systematic screen such as PACT *premier*[™]. The two screens provide a thorough exploration of crystallization conditions and the unique design of PACT *premier*[™] facilitates rational interpretation of results from both itself and JCSG-*plus*[™] assisting the design of subsequent experiments.

Formulation Notes:

JCSG-*plus* reagents are formulated using ultrapure water (>18.0 MΩ) and are sterile-filtered using 0.22 μm filters. No preservatives are added.

50% Stock solutions of Jeffamine are adjusted to pH 7.0 using HCl prior to inclusion in the reagents. Final pH may vary from that specified on the datasheet. Molecular Dimensions will be happy to discuss the precise formulation of individual reagents.

Individual reagents and stock solutions for optimization are available from Molecular Dimensions.

Enquiries regarding JCSG-*plus* formulation, interpretation of results or optimization strategies are welcome. Please e-mail, fax or phone your query to Molecular Dimensions.



Contact and product details can be found at www.moleculardimensions.com

References.

1. Page *et al* (2003). Shotgun crystallization strategy for structural genomics: an optimized two-tiered crystallization screen against the *Thermotoga maritima* proteome. *Acta Cryst.* **D59**, 1028-1037
2. Newman *et al* (2005). Towards rationalization of crystallization screening for small- to medium-sized academic laboratories: the PACT/JCSG+ strategy. *Acta Cryst.* **D61**, 1426-1431
3. McPherson *et al* (2001). A comparison of salts for the crystallisation of macromolecules, *Protein Science* **10**, 418422
4. Crystallization of Nucleic Acids and Proteins, Edited by A. Ducruix and R. Giegé, The Practical Approach Series, Oxford Univ. Press, 1992
5. Protein Crystallization Techniques Strategies & Tips, Edited by Terese Bergfors, IUL 1999
6. Methods and Results in the Crystallization of Membrane Proteins, Edited by So Iwata, IUL 2003.

Hints & Tips:

The JCSG-*plus*[™] sparse matrix screen is highly effective when used alongside a systematic screen such as PACT *premier*[™]. The two screens provide a thorough exploration of crystallization conditions and the unique design of PACT *premier*[™] facilitates rational interpretation of results from both itself and JCSG-*plus*[™] assisting the design of subsequent experiments.





JCSG-plus™ HT-96 Eco Screen

Wells A1 - D12

MD1-40-ECO

Well #	Conc.	Salt	Conc.	Buffer	pH	Conc.	Precipitant
A1	0.2 M	Lithium sulfate	0.1 M	Sodium acetate	4.5	50 % w/v	PEG 400
A2		None	0.1 M	Sodium citrate	5.5	20 % w/v	PEG 3000
A3	0.2 M	Ammonium citrate dibasic		None	-	20 % w/v	PEG 3350
A4	0.02 M	Calcium chloride dihydrate	0.1 M	Sodium acetate	4.6	30 % v/v	MPD
A5	0.2 M	Magnesium formate dihydrate		None	-	20 % w/v	PEG 3350
A6	0.2 M	Lithium sulfate	0.1 M	Phosphate/citrate	4.2	20 % w/v	PEG 1000
A7		None	0.1 M	CHES	9.5	20 % w/v	PEG 8000
A8	0.2 M	Ammonium formate		None	-	20 % w/v	PEG 3350
A9	0.2 M	Ammonium chloride		None	-	20 % w/v	PEG 3350
A10	0.2 M	Potassium formate		None	-	20 % w/v	PEG 3350
A11	0.2 M	Ammonium phosphate monobasic	0.1 M	Tris	8.5	50 % v/v	MPD
A12	0.2 M	Potassium nitrate		None	-	20 % w/v	PEG 3350
B1	0.8 M	Ammonium sulfate	0.1 M	Citrate	4.0		None
B2	0.2 M	Sodium thiocyanate		None	-	20 % w/v	PEG 3350
B3		None	0.1 M	BICINE	9.0	20 % w/v	PEG 6000
B4		None	0.1 M	HEPES	7.5	10 % w/v	PEG 8000
B5		None	0.1 M	MES	6.5	40 % v/v	MPD
B6		None	0.1 M	Phosphate/citrate	4.2	40 % v/v	Ethanol
B7		None	0.1 M	Sodium acetate	4.6	8 % w/v	PEG 4000
B8	0.2 M	Magnesium chloride hexahydrate	0.1 M	Tris	7.0	10 % w/v	PEG 8000
B9		None	0.1 M	Citrate	5.0	20 % w/v	PEG 6000
B10	0.2 M	Magnesium chloride hexahydrate	0.1 M	MES	6.5	50 % v/v	PEG 200
B11	1.6 M	Sodium citrate tribasic dihydrate		None	6.5		None
B12	0.2 M	Potassium citrate tribasic monohydrate		None	-	20 % w/v	PEG 3350
C1	0.2 M	Sodium chloride	0.1 M	Phosphate/citrate	4.2	20 % w/v	PEG 8000
C2	1.0 M	Lithium chloride	0.1 M	Citrate	4.0	20 % w/v	PEG 6000
C3	0.2 M	Ammonium nitrate		None	-	20 % w/v	PEG 3350
C4		None	0.1 M	HEPES	7.0	10 % w/v	PEG 6000
C5	0.8 M	Sodium phosphate monobasic monohydrate	0.1 M	Sodium HEPES	7.5		None
C6		None	0.1 M	Phosphate/citrate	4.2	40 % v/v	PEG 300
C7	0.2 M	Zinc acetate dihydrate	0.1 M	Sodium acetate	4.5	10 % w/v	PEG 3000
C8		None	0.1 M	Tris	8.5	20 % v/v	Ethanol
C9		None	0.1 M	Sodium/potassium phosphate	6.2	25 % v/v	1,2-Propandiol
C10		None	0.1 M	BICINE	9.0	10 % w/v	PEG 20,000
C11	2.0 M	Ammonium sulfate	0.1 M	Sodium acetate	4.6		None
C12		None		None	-	10 % w/v	PEG 1000
D1		None		None	-	24 % w/v	PEG 1500
D2	0.2 M	Magnesium chloride hexahydrate	0.1 M	Sodium HEPES	7.5	30 % v/v	PEG 400
D3	0.2 M	Sodium chloride	0.1 M	Sodium/potassium phosphate	6.2	50 % v/v	PEG 200
D4	0.2 M	Lithium sulfate	0.1 M	Sodium acetate	4.5	30 % w/v	PEG 8000
D5		None	0.1 M	HEPES	7.5	70 % v/v	MPD
D6	0.2 M	Magnesium chloride hexahydrate	0.1 M	Tris	8.5	20 % w/v	PEG 8000
D7	0.2 M	Lithium sulfate	0.1 M	Tris	8.5	40 % v/v	PEG 400
D8		None	0.1 M	Tris	8.0	40 % v/v	MPD
D9	0.17 M	Ammonium sulfate		None	-	26 % w/v	PEG 4000
D10		None		None	-	15 % v/v	Glycerol
D10	0.2 M	Calcium acetate hydrate	0.1 M	MES	6.5	40 % v/v	PEG 300
D11	0.14 M	Calcium chloride dihydrate	0.1 M	Sodium acetate	4.6	14 % v/v	2-Propanol
D11		None		None	-	30 % v/v	Glycerol
D12	0.04 M	Potassium phosphate monobasic		None	-	16 % w/v	PEG 8000
D12		None		None	-	20 % v/v	Glycerol





JCSG-plus™ HT-96 Eco Screen Wells E1 - H12

MD1-40-ECO

Well #	Conc.	Salt	Conc.	Buffer	pH	Conc.	Precipitant
E1	1.0 M	Sodium citrate tribasic dihydrate	0.1 M	MES	6.5		None
E2	2.0 M	Ammonium sulfate	0.1 M	MES	6.5		None
	0.2 M	Sodium chloride					
E3	0.2 M	Sodium chloride	0.1 M	HEPES	7.5	10 % v/v	2-Propanol
E4	1.26 M	Ammonium sulfate	0.1 M	Tris	8.5		None
	0.2 M	Lithium sulfate					
E5		None	0.1 M	CAPS	10.5	40 % v/v	MPD
E6	0.2 M	Zinc acetate dihydrate	0.1 M	Imidazole	8.0	20 % w/v	PEG 3000
E7	0.2 M	Zinc acetate dihydrate	0.1 M	MES	6.5	10 % v/v	2-Propanol
E8	1.0 M	Ammonium phosphate dibasic	0.1 M	Sodium acetate	4.5		None
E9	1.6 M	Magnesium sulfate heptahydrate	0.1 M	MES	6.5		None
E10		None	0.1 M	BICINE	9.0	10 % w/v	PEG 6000
E11	0.16 M	Calcium acetate hydrate	0.08 M	MES	6.5	14.4 % w/v	PEG 8000
						20 % v/v	Glycerol
E12		None	0.1 M	Imidazole	8.0	10 % w/v	PEG 8000
F1	0.05 M	Cesium chloride	0.1 M	MES	6.5	30 % v/v	Jeffamine® M-600
F2	3.2 M	Ammonium sulfate	0.1 M	Citrate	5.0		None
F3		None	0.1 M	Tris	8.0	20 % v/v	MPD
F4		None	0.1 M	HEPES	7.5	20 % v/v	Jeffamine® M-600
F5	0.2 M	Magnesium chloride hexahydrate	0.1 M	Tris	8.5	50 % v/v	Ethylene glycol
F6		None	0.1 M	BICINE	9.0	10 % v/v	MPD
F7	0.8 M	Succinic acid pH 7.0		None			None
F8	2.1 M	DL-Malic acid pH 7.0		None			None
F9	2.4 M	Sodium malonate dibasic monohydrate pH 7.0		None			None
F10	1.1 M	Sodium malonate dibasic monohydrate	0.1 M	HEPES	7.0	0.5 % v/v	Jeffamine® ED-2003
F11	1.0 M	Succinic acid	0.1 M	HEPES	7.0	1 % w/v	PEG 2000 MME
F12		None	0.1 M	HEPES	7.0	30 % v/v	Jeffamine® M-600
G1		None	0.1 M	HEPES	7.0	30 % v/v	Jeffamine® ED-2003
G2	0.02 M	Magnesium chloride hexahydrate	0.1 M	HEPES	7.5	22 % w/v	Poly(acrylic acid sodium salt) 5100
G3	0.01 M	Cobalt(II) chloride hexahydrate	0.1 M	Tris	8.5	20 % w/v	Polyvinylpyrrolidone
G4	0.2 M	TMAO	0.1 M	Tris	8.5	20 % w/v	PEG 2000 MME
G5	0.005 M	Cobalt(II) chloride hexahydrate	0.1 M	HEPES	7.5	12 % w/v	PEG 3350
	0.005 M	Cadmium chloride hemi(pentahydrate)					
	0.005 M	Magnesium chloride hexahydrate					
	0.005 M	Nickel(II) chloride hexahydrate					
G6	0.2 M	Sodium malonate dibasic monohydrate		None		20 % w/v	PEG 3350
G7	0.1 M	Succinic acid		None		15 % w/v	PEG 3350
G8	0.15 M	DL-Malic acid		None		20 % w/v	PEG 3350
G9	0.1 M	Potassium thiocyanate		None		30 % w/v	PEG 2000 MME
G10	0.15 M	Potassium bromide		None		30 % w/v	PEG 2000 MME
G11	2.0 M	Ammonium sulfate	0.1 M	BIS-Tris	5.5		None
G12	3.0 M	Sodium chloride	0.1 M	BIS-Tris	5.5		None
H1	0.3 M	Magnesium formate dihydrate	0.1 M	BIS-Tris	5.5		None
H2	1.0 M	Ammonium sulfate	0.1 M	BIS-Tris	5.5	1 % w/v	PEG 3350
H3		None	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H4	0.2 M	Calcium chloride dihydrate	0.1 M	BIS-Tris	5.5	45 % v/v	MPD
H5	0.2 M	Ammonium acetate	0.1 M	BIS-Tris	5.5	45 % v/v	MPD
H6	0.1 M	Ammonium acetate	0.1 M	BIS-Tris	5.5	17 % w/v	PEG 10,000
H7	0.2 M	Ammonium sulfate	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H8	0.2 M	Sodium chloride	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H9	0.2 M	Lithium sulfate	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H10	0.2 M	Ammonium acetate	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H11	0.2 M	Magnesium chloride hexahydrate	0.1 M	BIS-Tris	5.5	25 % w/v	PEG 3350
H12	0.2 M	Ammonium acetate	0.1 M	HEPES	7.5	45 % v/v	MPD

Abbreviations: Bis-Tris; Bis-(2-hydroxyethyl)imino-tris(hydroxymethyl)methane, **CAPS**; N-Cyclohexyl-3-aminopropanesulfonic acid, **CHES**; 2-(N-Cyclohexylamino)ethane Sulfonic Acid, **HEPES**; 2-(4-(2-Hydroxyethyl)-1-piperazinyl)ethanesulfonic Acid, **Sodium HEPES**; 2-(4-(2-Hydroxyethyl)-1-piperazinyl)ethanesulfonic Acid Sodium Salt, **MES**; 2-(N-morpholino)ethanesulfonic acid, **MPD**; 2,4-methyl pentanediol, **PEG**; Polyethylene glycol, **TMAO**: Trimethylamine N-oxide, **Tris**; 2-Amino-2-(hydroxymethyl)propane-1,3-diol.

N.B. Jeffamine® ED-2001 has been superseded with Jeffamine® ED-2003. Polyvinylpyrrolidone K15 is called Polyvinylpyrrolidine.



Manufacturer's safety data sheets are available from our website or by scanning the QR code here:



Re-Ordering details:

Catalogue Description Code	Pack size	Catalogue
JCSG- <i>plus</i> TM	96 x 10 mL	MD1-37
JCSG- <i>plus</i> TM HT-96	96 x 1 mL	MD1-40
Eco Screens		
JCSG- <i>plus</i> TM Eco Screen	96 x 10 mL	MD1-37-ECO
JCSG- <i>plus</i> TM HT-96 Eco Screen	96 x 1 mL	MD1-40-ECO
Green Screens (contain fluorescent green dye - ideal for UV)		
JCSG- <i>plus</i> TM Green Screen	96 x 10 mL	MD1-56
JCSG- <i>plus</i> TM HT-96 Green Screen	96 x 1 mL	MD1-53
Combo Packs		
Super2 Combo Value Pack (JCSG- <i>plus</i> TM + PACT <i>premier</i> TM)	2 x 96 x 10 mL	MD1-75
Super2 Combo HT-96 Value Pack (JCSG- <i>plus</i> TM HT-96 + PACT <i>premier</i> TM HT-96)	2 x 96 x 10 mL	MD1-75-HT
Single Reagents		
JCSG- <i>plus</i> TM single reagents	100 mL	MDSR-37-tube number
JCSG- <i>plus</i> TM single reagents	100 mL	MDSR-40-well number

For JCSG-*plus*TM stock solutions please visit the Optimization section on our website.